

Home // Tracheal Laryngeal Tubes // ENT // JACKSON IMPROVED STAINLESS STEEL TRACHEA TUBE #5 FOR SHIKANI-FRENCH SPEAKING VALVE

Catalog # :

518034A

JACKSON IMPROVED STAINLESS STEEL TRACHEA TUBE #5 FOR SHIKANI-FRENCH SPEAKING VALVE

Description:

6.2mm I.D., 9.0mm O.D., 68mm, fenestrated, inner and outer

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699
