

Home // Tracheal Laryngeal Tubes // ENT // JACKSON IMPROVED STAINLESS STEEL TRACHEA TUBE #6 FOR SHIKANI-FRENCH SPEAKING VALVE


Catalog # :

518039A

JACKSON IMPROVED STAINLESS STEEL TRACHEA TUBE #6 FOR SHIKANI-FRENCH SPEAKING VALVE

Description:

7.2mm I.D., 10.0mm O.D., 69mm, fenestrated, inner and outer



NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699
