

Home // Tracheal Laryngeal Tubes // ENT // SHIKANI-FRENCH SPEAKING VALVE FOR USE WITH JACKSON IMPROVED STAINLESS STEEL TRACHEOSTOMY TUBES

Catalog # :


**518085**

---

## SHIKANI-FRENCH SPEAKING VALVE FOR USE WITH JACKSON IMPROVED STAINLESS STEEL TRACHEOSTOMY TUBES

**Description:**

package of 2



---

United States

To order this product, please call:

**866-246-6990**

---

Canada

To order this product, please call:

**800-387-9699**

---