

Home // Tracheal Laryngeal Tubes // ENT // JACKSON IMPROVED TRACHEOSTOMY TUBE #5, W/15MM ADAPTER, REGULAR LENGTH

Catalog # :

~~518431~~

JACKSON IMPROVED TRACHEOSTOMY TUBE #5, W/15MM ADAPTER, REGULAR LENGTH

Description:

6.2mm I.D., 9mm O.D., 68mm, stainless steel

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699