

Home // Tracheal Laryngeal Tubes // ENT // JACKSON ORIGINAL TRACHEOSTOMY TUBE #4, W/15MM ADAPTER, REGULAR LENGTH

Catalog # :

**518526**

---

## JACKSON ORIGINAL TRACHEOSTOMY TUBE #4, W/15MM ADAPTER, REGULAR LENGTH

Description:

5.5mm I.D., 8mm O.D., 62mm, stainless steel

NO IMAGE  
AVAILABLE

---

United States

To order this product, please call:

**866-246-6990**

---

Canada

To order this product, please call:

**800-387-9699**

---