

Home // Tracheal Laryngeal Tubes // ENT // JACKSON ORIGINAL TRACHEOSTOMY TUBE #6, W/15MM ADAPTER, REGULAR LENGTH

Catalog # :

~~518536~~

JACKSON ORIGINAL TRACHEOSTOMY TUBE #6, W/15MM ADAPTER, REGULAR LENGTH

Description:

7.2mm I.D., 10mm O.D., 69mm, stainless steel

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699