

Home // Tracheal Laryngeal Tubes // ENT // JACKSON ORIGINAL TRACHEOSTOMY TUBE #7, W/15MM ADAPTER, REGULAR LENGTH

Catalog # :

~~518541~~

JACKSON ORIGINAL TRACHEOSTOMY TUBE #7, W/15MM ADAPTER, REGULAR LENGTH

Description:

8.1mm I.D., 11mm O.D., 69mm, stainless steel

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699
