

Home // Tracheal Laryngeal Tubes // ENT // JACKSON IMPROVED TRACHEOSTOMY TUBE #5, W/15MM ADAPTER, SHORT LENGTH

Catalog # :

~~518625~~

---

## JACKSON IMPROVED TRACHEOSTOMY TUBE #5, W/15MM ADAPTER, SHORT LENGTH

Description:

6.2mm I.D., 9mm O.D., 56mm, stainless steel



---

United States

To order this product, please call:

**866-246-6990**

---

Canada

To order this product, please call:

**800-387-9699**

---