

Home // Tracheal Laryngeal Tubes // ENT // JACKSON IMPROVED TRACHEOSTOMY TUBE #6, W/15MM ADAPTER, SHORT LENGTH

Catalog # :

~~518630~~

JACKSON IMPROVED TRACHEOSTOMY TUBE #6, W/15MM ADAPTER, SHORT LENGTH

Description:

7.2mm I.D., 10mm O.D., 56mm, stainless steel



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699
