

Home // Tracheal Laryngeal Tubes // ENT // JACKSON IMPROVED TRACHEOSTOMY TUBE #7, W/15MM ADAPTER, SHORT LENGTH

Catalog # :

~~518635~~

JACKSON IMPROVED TRACHEOSTOMY TUBE #7, W/15MM ADAPTER, SHORT LENGTH

Description:

8.1mm I.D., 11mm O.D., 56mm, stainless steel

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699
