

Home // Tracheal Laryngeal Tubes // ENT // JACKSON IMPROVED TRACHEOSTOMY TUBE #8, W/15MM ADAPTER, SHORT LENGTH

Catalog # :

~~518640~~

## JACKSON IMPROVED TRACHEOSTOMY TUBE #8, W/15MM ADAPTER, SHORT LENGTH

Description:

9.2mm I.D., 12mm O.D., 56mm, stainless steel

NO IMAGE  
AVAILABLE

United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**