

Home // Tracheal Laryngeal Tubes // ENT // JACKSON ORIGINAL TRACHEOSTOMY TUBE #4, W/15MM ADAPTER, SHORT LENGTH

Catalog # :
~~518720~~

JACKSON ORIGINAL TRACHEOSTOMY TUBE #4, W/15MM ADAPTER, SHORT LENGTH

Description:

5.5mm I.D., 8mm O.D., 52mm, stainless steel

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699