

Home // Laryngoscopes, Bronchoscopes, Esophagoscopes, Fiberoptic Cables, Light Carriers // Laryngoscopes // ENT
// HOLINGER LARYNGOSCOPE, INFANT

Catalog # :

~~521500~~

HOLINGER LARYNGOSCOPE, INFANT

Description:

Overall length 82mm, min. I.D. 8.0mm, replacement light carrier 526204

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699