

Home // Laryngoscopes, Bronchoscopes, Esophagoscopes, Fiberoptic Cables, Light Carriers // Laryngoscopes // ENT
// HOLINGER ANTERIOR COMMISSURE LARYNGOSCOPE, ADULT


Catalog # :

522010

HOLINGER ANTERIOR COMMISSURE LARYNGOSCOPE, ADULT

Description:

Overall length 170mm, min. I.D. 8.3mm, replacement light carrier 526028



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699
