

Home // Laryngoscopes, Bronchoscopes, Esophagoscopes, Fiberoptic Cables, Light Carriers // Laryngoscopes // ENT  
// HOLINGER ANTERIOR COMMISSURE LARYNGOSCOPE, ADULT EXTRA-LONG

Catalog # :

~~522012~~

---

## HOLINGER ANTERIOR COMMISSURE LARYNGOSCOPE, ADULT EXTRA-LONG

**Description:**

Overall length 180mm, min. I.D. 8.3mm, replacement light carrier 526038



---

United States

To order this product, please call:

**866-246-6990**

---

Canada

To order this product, please call:

**800-387-9699**

---