

Home // Laryngoscopes, Bronchoscopes, Esophagoscopes, Fiberoptic Cables, Light Carriers // Laryngoscopes // ENT
// HOLINGER "HOUR GLASS" ANTERIOR COMMISSURE LARYNGOSCOPE, ADULT EXTRA LONG

Catalog # :

522032

HOLINGER "HOUR GLASS" ANTERIOR COMMISSURE LARYNGOSCOPE, ADULT EXTRA LONG

Description:

Overall length 181mm, min. I.D. 7.3mm, replacement light carrier 526214



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699
