

Catalog # :

522190**Ossoff-Pilling® Micro-Laryngoscope, Adult Female****Description:**


The new Ossoff-Pilling Micro-Laryngoscope is available in two proximal widths. The 25.5 mm width (522191) can be introduced in most patients and is easily aligned with the microscope (400 mm objective). The 19 mm scope (522190) has the narrowest width capable of absorbing both microscope sight lines and can be used in most patients who previously could accept only the original Holinger scope.

- Holinger distal end, ease of introduction, excellent visualization.
- 19 or 25.5 mm proximal widths, increases patient selection for micro procedures.
- Standard adult handle, can be used with most suspension systems, B.U. suspension handle available as option.
- Single light carrier, external tube dimensions are kept to a minimum.
- Laser finish, laser safety.
- Supplied with one light carrier. Cable sold separately

522190: Ossoff-Pilling Micro-Laryngoscope, Adult Female Overall length **168 mm**, Min Distal Opening **7.0 mm**, Nominal Proximal Dims **19.0 mm x 15.0 mm**, Replace Carrier 527864

Related:

522191: 24mm wide distal opening



United States

Details:

To order this product, please call:

- 19mm wide distal opening

866-246-6990

Canada

To order this product, please call:

800-387-9699