

Catalog # :

522221**Low-Profile Dedo-Style Laser Laryngoscope, Adult****Description:**

The Newest Addition to the Pilling® Family. Featuring a low-profile design to accommodate a wider range of patient anatomies than previous-generation Pilling Dedo laryngoscopes, while maintaining the same internal diameter for binocular visualization. The Low-Profile Dedo-Style Laser Laryngoscope has an outer diameter that is 20% (6.5mm) slimmer than our Dedo Laser Laryngoscope 522228, and 5% (1.5mm) slimmer than our Dedo Micro Laryngoscope 522225P. Designed to be safer when using laser, this new laryngoscope contains satin finish outside, matte finish inside and a smoke evacuation channel.


522221: Low-Profile Dedo-Style Laser Laryngoscope, Adult is supplied with a single light carrier, flush with the handle, which accommodates a slim and lightweight blue Pilling fiberoptic cable (sold separately). Size **Adult**. Length **178mm**, Min I.D. **10.7mm**, Nominal Proximal Dims **19.0mm x 19.0mm**, Max O.D. **25.5mm**, Replace Carrier 527868

Related:

522225P: Dedo-Pilling Micro-Laryngoscope, Adult

522228: Dedo Laser Micro Laryngoscope, Adult

522231: Ossoff-Karlan-Dedo Laser Laryngoscope, Adult



To order this product, please call:

+1 (919) 433-4999