


Catalog # :

**522300P****CARLENS MEDIASTINOSCOPE, SHORT****Description:**

The Carlens Mediastinoscopes come in short, long and extra-long lengths. The short model is suitable for inspection of the mediastinum above the bifurcation, while the long and extra-long models are designed for the region of the carina. Supplied with one light carrier.

Short scope, overall length 116 mm, min. I.D. 10.8 mm, nominal proximal dims 16.0 mm, max O.D. 18.4 mm, replacement light carrier 526350



United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**