


Catalog # :

~~522302~~**CARLENS MEDIASTINOSCOPE, EXTRA-LONG****Description:**

The Carlens Mediastinoscopes come in short, long and extra-long lengths. The short model is suitable for inspection of the mediastinum above the bifurcation, while the long and extra-long models are designed for the region of the carina. Supplied with one light carrier.

Extra-long scope, overall length 191 mm, min. I.D. 11.9 mm, nominal proximal dims 19.0 mm, max O.D. 22.3 mm, replacement light carrier 526352



United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**