

Home // Laryngoscopes, Bronchoscopes, Esophagoscopes, Fiberoptic Cables, Light Carriers // Bronchoscopes // ENT
// HOLINGER BRONCHOSCOPE 3.5X30

Catalog # :

~~522577~~

HOLINGER BRONCHOSCOPE 3.5X30

Description:

Scope size 3.5 x 30, overall length 305mm, min. I.D. 3.4mm, replacement light carrier 527510

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699