

Home // Laryngoscopes, Bronchoscopes, Esophagoscopes, Fiberoptic Cables, Light Carriers // Bronchoscopes // ENT  
// HOLINGER BRONCHOSCOPE 5X30

Catalog # :

~~522579~~

## HOLINGER BRONCHOSCOPE 5X30

### Description:

Scope size 5 x 30, overall length 305mm, min. I.D. 4.5mm, replacement light carrier 527510

NO IMAGE  
AVAILABLE

United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**