

Home // Laryngoscopes, Bronchoscopes, Esophagoscopes, Fiberoptic Cables, Light Carriers // Bronchoscopes // ENT
// HOLINGER BRONCHOSCOPE 5X45

Catalog # :

~~522582~~

HOLINGER BRONCHOSCOPE 5X45

Description:

Scope size 5 x 45, overall length 451mm, min. I.D. 4.9mm, replacement light carrier 527525

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699