

Home // Laryngoscopes, Bronchoscopes, Esophagoscopes, Fiberoptic Cables, Light Carriers // Bronchoscopes // ENT
// HOLINGER BRONCHOSCOPE 9X40

Catalog # :

~~522587~~

HOLINGER BRONCHOSCOPE 9X40

Description:

Scope size 9 x 40, overall length 400mm, min. I.D. 8.5mm, replacement light carrier 527580

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699