

Home // Laryngoscopes, Bronchoscopes, Esophagoscopes, Fiberoptic Cables, Light Carriers // Esophagoscopes // ENT
// HOLINGER ESOPHAGOSCOPE

Catalog # :

~~523582~~

HOLINGER ESOPHAGOSCOPE

Description:

9 x 45, Overall length 451mm, min. I.D. 8.5mm, replacement light carrier 527583

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699