

Home // Laryngoscopes, Bronchoscopes, Esophagoscopes, Fiberoptic Cables, Light Carriers // Esophagoscopes // ENT
// JACKSON FIBEROPTIC OVAL ESOPHAGOSCOPE

Catalog # :

~~523600P~~

JACKSON FIBEROPTIC OVAL ESOPHAGOSCOPE

Description:

7 x 10 x 25, Overall length 251mm, min. I.D. 6.4mm, nominal proximal dims 7.0 X 10.0mm, replacement light carrier 526262

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699