

Catalog # :**806034****Multipurpose Stainless Steel Tray with Lid**

This item has been discontinued and is no longer available.

Description:

Removable locking lid stores easily under tray during use.

Stainless steel storage tray with lid, 3-1/4" x 7" x 1-3/4" (8.5 x 18 x 4.5 cm)

Related:

806040 Multipurpose Stainless Steel Tray with Lid, 6" x 30" x 3-1/2"

806048 Multipurpose Stainless Steel Tray with Lid, 15" x 10-1/2" x 3-1/2"

806050 Multipurpose Stainless Steel Tray with Lid, 20" x 10-1/2" x 3-1/2"

806052 Multipurpose Stainless Steel Tray with Lid, 20-1/4" x 13" x 3-1/2"



To order this product, please call:

+1 (919) 433-4999

Details:

- Removable locking lid stores easily under tray during use.
- Stainless steel storage tray with lid, 3-1/4" x 7" x 1-3/4" (8.5 x 18 x 4.5 cm)