

Home // Bending Irons, Rod Holders, Ligature Carriers, Wire Passers // Bending Irons, Rod Holders // Orthopedic // WIRE AND PIN TRACTOR (BOEHLER) - LARGE W/4.0mm STEIN PIN

Catalog # :

~~KM48330~~

## WIRE AND PIN TRACTOR (BOEHLER) - LARGE W/4.0mm STEIN PIN

### Description:

Large, 210 mm inside length, 150 mm inside width. Pin size and diameter: 9" 4.0 mm.  
Includes pin KM71525. 6/pack.

NO IMAGE  
AVAILABLE

United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**