

Catalog # :


~~KM98150~~**HOKE SMARTRACK OSTEOTOME****Description:**

- Nylon coated brackets help avoid metal-to-metal contact.
- Plastic plate on underside of top bracket helps protect sharp tips from being damaged.
- Side arm stays locked to help prevent instruments from falling out, even when turned upside down.
- Deep brackets hold instruments in place.
- Custom designed base keeps instruments upright during loading and unloading.
- Reinforced, sturdier frames are made from corrosion resistant anodized aluminum.
- Fully autoclavable to 270°F/132°C.
- Quick reference chart easily identifies the type and size instrument rack holds.
- Manufactured only by Teleflex®.

8-1/4" x 5" x 2-1/4"

(21 cm x 12.5 cm x 5.5 cm)

Holds eight 5-1/4" (13.5 cm) KMedic® Hoke Chisels or Osteotomes



United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**