

Catalog # :

W356PXA

Pilling® Bifurcated Endfitting Blue Fiberoptic Cable

Description:

For use with Pilling® Laryngoscopes, Bronchoscopes and Esophagoscopes.

W356PXA: Pilling® Blue Fiber Optic Cable , Bifurcated cable, Wolf illuminator end, 3.5 mm bundle, 6 ft. length, Pilling® rotating instrument end

Related:

A358AFA: ACMI illuminator end, 8 ft. length, ACMI female instrument end

A358PRA: ACMI illuminator end, 8 ft. length, Pilling® rotating instrument end

O358PRA: Olympus illuminator end, 8 ft. length, Pilling® rotating instrument end

P356PRA: Pilling® illuminator end, 6 ft. length, Pilling® rotating instrument end

P358PRA: Pilling® illuminator end, 8 ft. length, Pilling® rotating instrument end

S358PRA: Storz illuminator end, 8 ft. length, Pilling® rotating instrument end

W356PXA: Bifurcated cable, Wolf illuminator end, 6 ft. length, Pilling® rotating instrument end

W358PRA: Wolf illuminator end, 8 ft. length, Pilling® rotating instrument end

W358PXA: Bifurcated cable, Wolf illuminator end, 8 ft. length, Pilling® rotating instrument end



United States

To order this product, please call:

866-246-6990

Details:

- For use with Pilling® Laryngoscopes, Bronchoscopes and

Canada

To order this product, please call:

800-387-9699

Esophagoscopes.

Bifurcated cable, Wolf illuminator end, 3.5 mm bundle, 6 ft. length, Pilling® rotating instrument end