

Home // Mirrors, Snares, Mouth Gags

**Mirrors, Snares, Mouth Gags****461126 : MCIVOR FRAME****508200 : LARYNGEAL MIRROR**

#000, 8mm diameter, Shaft length is 6 3/4" (17.0cm)

**508201 : LARYNGEAL MIRROR**


#00, 10mm diameter, Shaft length is 6 3/4" (17.0cm)

**508202 : LARYNGEAL MIRROR**

#0, 12mm diameter, Shaft length is 6 3/4" (17.0cm)

**508203 : LARYNGEAL MIRROR**

#1, 16mm diameter, Shaft length is 6 3/4" (17.0cm)

**508204 : LARYNGEAL MIRROR**

#2, 18mm diameter, Shaft length is 6 3/4" (17.0cm)



**508205 : LARYNGEAL MIRROR**

#3, 20mm diameter, Shaft length is 6 3/4" (17.0cm)



**508206 : LARYNGEAL MIRROR**

#4, 22mm diameter, Shaft length is 6 3/4" (17.0cm)



**508207 : LARYNGEAL MIRROR**

#5, 24mm diameter, Shaft length is 6 3/4" (17.0cm)




**508208 : LARYNGEAL MIRROR**

#6, 28mm diameter, Shaft length is 6 3/4" (17.0cm)



**508209 : LARYNGEAL MIRROR**

#7, 30mm diameter, Shaft length is 6 3/4" (17.0cm)



**508230 : JAKO MIRROR**

Working length 8 1/4" (21.0cm)



**509209 : LASER MIRROR**

Highly polished stainless front surface, 6mm diameter, Working length 8 1/4" (21.0cm)

**509302 : FRONT SURFACE MIRROR**

3mm diameter, Working length 11" (28.0cm)



**P24263 : WILDER LOOP LENS**

Serrated

NO IMAGE  
AVAILABLE

**P24270 : KIRBY EXPRESSOR HOOK AND LOOP**

NO IMAGE  
AVAILABLE

**P24279 : FOSTER ENUCLEATION SNARE**

w/ pkg. of 12 snares

NO IMAGE  
AVAILABLE

**P24280 : TWISTED SNARE WIRE**

pkg. of 12

NO IMAGE  
AVAILABLE

**P26434 : MIDDLE EAR MIRROR**

Angled 45 deg., 2mm, 6 1/2" (16.5cm)

NO IMAGE  
AVAILABLE

**P28275 : OCTAGONAL CONE SOCKET MIRROR HANDLE**

4 1/2" (11.5cm)

NO IMAGE  
AVAILABLE

AVAILABLE

NO IMAGE  
AVAILABLE

**P28319 : SLUDER JANSEN MOUTH GAG, CHILD**

NO IMAGE  
AVAILABLE

**P28321 : SLUDER JANSEN MOUTH GAG, ADULT**

NO IMAGE  
AVAILABLE

**P28325 : MCIVOR MOUTH GAG, TONGUE BLADES, SIZE #1**

22mm x 69mm

NO IMAGE  
AVAILABLE

**P28326 : MCIVOR MOUTH GAG, TONGUE BLADES, SIZE #2**

24mm x 90mm

NO IMAGE  
AVAILABLE

**P28327 : MCIVOR MOUTH GAG, TONGUE BLADES, SIZE #3**

25mm x 100mm