

Laryngoscopes, Bronchoscopes, Esophagoscopes, Fiberoptic Cables, Light Carriers**522584 : HOLINGER BRONCHOSCOPE 6X40**

Scope size 6 x 40, overall length 400mm, min. I.D. 5.6mm, replacement light carrier 527560

**522585 : HOLINGER BRONCHOSCOPE 7X40**

Scope size 7 x 40, overall length 400mm, min. I.D. 5.8mm, replacement light carrier 527560

**522586 : HOLINGER BRONCHOSCOPE 8X40**

Scope size 8 x 40, overall length 400mm, min. I.D. 7.6mm, replacement light carrier 527580

**522587 : HOLINGER BRONCHOSCOPE 9X40**

Scope size 9 x 40, overall length 400mm, min. I.D. 8.5mm, replacement light carrier 527580

**523120P : LASER BRONCHOSCOPE 5X30**

Scope size 5.0 x 30, overall length 300mm, min. I.D. 5.1mm, replacement light carrier 527550

**523130P : LASER BRONCHOSCOPE 6X35**

Scope size 6.0 x 35, overall length 350mm, min. I.D. 6.0mm, replacement light carrier 523131



523131 : LIGHT CARRIER OVAL LARGE 35CM

For scope number 523130P



523140P : LASER BRONCHOSCOPE 7X40

Scope size 7.0 x 40, overall length 400mm, min. I.D. 7.0mm, replacement light carrier 527580



523150P : LASER BRONCHOSCOPE 8X40

Scope size 8.0 x 40, overall length 400mm, min. I.D. 8.0mm, replacement light carrier 527580



523160P : LASER BRONCHOSCOPE 9X40

Scope size 9.0 x 40, overall length 400mm, min. I.D. 9.0mm, replacement light carrier 527580



523505 : HOLINGER ESOPHAGOSCOPE

Infant, 6 x 20, Overall length 200mm, min. I.D. 4.0mm, replacement light carrier 526218



523570 : HOLINGER ESOPHAGOSCOPE

Infant, 3.5 x 30, Overall length 300mm, min. I.D. 3.4mm, replacement light carrier 527510



523571 : HOLINGER ESOPHAGOSCOPE

4 x 30, Overall length 300mm, min. I.D. 3.9mm, replacement light carrier 527510

523572 : HOLINGER ESOPHAGOSCOPE

6 x 30, Overall length 300mm, min. I.D. 5.4mm, replacement light carrier 527550



523573 : HOLINGER ESOPHAGOSCOPE

7 x 30, Overall length 300mm, min. I.D. 6.8mm, replacement light carrier 527550



523574 : HOLINGER ESOPHAGOSCOPE

9 x 30, Overall length 300mm, min. I.D. 8.8mm, replacement light carrier 527575



523575 : HOLINGER ESOPHAGOSCOPE

5 x 35, Overall length 349mm, min. I.D. 4.7mm, replacement light carrier 527515



523577 : HOLINGER ESOPHAGOSCOPE

7 x 45, Overall length 451mm, min. I.D. 6.8mm, replacement light carrier 527565



523579 : HOLINGER ESOPHAGOSCOPE


9 x 53, Overall length 530mm, min. I.D. 8.5mm, replacement light carrier 527585



523581 : HOLINGER ESOPHAGOSCOPE

8 x 45, Overall length 451mm, min. I.D. 7.6mm, replacement light carrier 527583





523582 : HOLINGER ESOPHAGOSCOPE

9 x 45, Overall length 451mm, min. I.D. 8.5mm, replacement light carrier 527583



523600P : JACKSON FIBEROPTIC OVAL ESOPHAGOSCOPE

7 x 10 x 25, Overall length 251mm, min. I.D. 6.4mm, nominal proximal dims 7.0 X 10.0mm, replacement light carrier 526262



523602 : JACKSON OVAL ESOPHAGOSCOPE

8 x 12 x 50, Overall length 502mm, min. I.D. 7.4mm, replacement light carrier 526284



523604 : JACKSON OVAL ESOPHAGOSCOPE

9 x 13 x 25, Overall length 251mm, min. I.D. 7.4mm, replacement light carrier 526262



523606 : JACKSON OVAL ESOPHAGOSCOPE

10 x 14 x 30, Overall length 302mm, min. I.D. 8.4mm, replacement light carrier 526266