

Laryngoscopes, Bronchoscopes, Esophagoscopes, Fiberoptic Cables, Light Carriers



526038 : LIGHT CARRIER TYPE 18H

For scope number 502022



526046 : LIGHT CARRIER TYPE 22

For scope number 523677



526049 : LIGHT CARRIER TYPE 25



526204 : LIGHT CARRIER TYPE 8PH

For scope number 521500



526206 : LIGHT CARRIER TYPE 9PH

For scope number 522002



526210 : LIGHT CARRIER TYPE 12PH

For scope number 522004



526212 : LIGHT CARRIER TYPE 17HG

For scope number 522045



526214 : LIGHT CARRIER TYPE 18HG

For scope number 522032



526218 : LIGHT CARRIER TYPE 20PH

For scope number 522500 / 523505



526262 : LIGHT CARRIER TYPE 25R

For scope number 523600 / 523604




526263 : LIGHT CARRIER TYPE 27R



526264 : LIGHT CARRIER TYPE 29J

For scope number 523658



526266 : LIGHT CARRIER TYPE 30R

For scope number 523606

526268 : LIGHT CARRIER TYPE 32J

For scope number 523650



526270 : LIGHT CARRIER TYPE 35J

For scope number 523652



526272 : LIGHT CARRIER TYPE 35R

For scope number 523614



526274 : LIGHT CARRIER TYPE 40R

For scope number 523608



526276 : LIGHT CARRIER TYPE 42J

For scope number 523654



526278 : LIGHT CARRIER TYPE 45J


For scope number 523656



526282 : LIGHT CARRIER TYPE 50J

For scope number 523660





526284 : LIGHT CARRIER TYPE 50R

For scope number 523602 / 523610 / 523612 / 523635



526350 : LIGHT CARRIER

For scope number 522300 / 522300P



526351 : LIGHT CARRIER

For scope number 522301



526352 : LIGHT CARRIER

For scope number 522302



526360 : LIGHT CARRIER LEFT

For scope number 522305 / 522305P

