

**Miscellaneous****KM49632 : HANDLE ONLY FOR KM 49630**

7 3/8" (18.5cm)

**89-2719 : SAPHLITE HANDLE, MODULAR****507485S : MICRO ATOMIZER (SM #'S 151 & 1-RD)****KM46892 : MURPHY BONE SKID**

11 3/4" (30.0cm)

**502574 : TYGON TUBING 4MM W/FITTING****KM47904 : SMALL QUICK COUPLING HANDLE**

for small and large screw sets, 4 1/2" (11.5cm)



**140879 : Modular Reusable Verres Body/Obturator 15 cm**



**507335 : BONNIER CATHETER ADAPTER**

1 3/8" (3.5cm)



**509216 : OPERATING PLATFORM**

Angled, Ebonized, 0.3" (8mm), 11 5/8" (29.5cm)




**509207 : VOCAL CORD ROLLERS**

Ebonized




**140880 : Modular Reusable Verres Body/Obturator 18 cm**



**P26511 : PICK**

3mm curve from tip, 6 5/16" (16cm)



**P30366 : CORKSCREW T-HANDLE**

7 1/2" (19cm)

**KM47902 : Mini Quick Coupling Handle**

For holding quick coupling drills, countersinks, taps and screwdrivers. For mini screw set. 4-1/8" (10.5 cm)



**89-2401 : PITIE-SALPETRIERE VEIN RTR HK 24CM, VEIN HOOK**

For use with surgical blades no. 10, 11, 12, and 15.



**411166 : FARRELL APPLICATOR TRIANGULAR**

6 3/8" (16.0cm)



**KM47906 : BALL SPIKE PELVIC STRAIGHT**

11 3/4" (30.0cm)



**KM46602 : PERCY AMPUTATION SHIELD**

12" (30.5cm)




**522590 : BRONCHOSCOPE ADAPTER 11MM TO 15MM**

3/4" (2.0cm)



**728167 : SPATULA ELECTRODE, 33CM**



AVAILABLE

NO IMAGE  
AVAILABLE

**P26601 : SHEEHY FASCIA PRESS**

23mm x 31mm plate, 7 1/2" (19cm)

NO IMAGE  
AVAILABLE

**728195 : APPEL TYPE KNOT PUSHER 5MM 40CM**

NO IMAGE  
AVAILABLE

**P17773 : VARADY EXTRACTOR**

Double-ended, 7" (18 cm). Crochet hook on one end.

NO IMAGE  
AVAILABLE

**P24338 : OLK RETINAL PICK**

Angled 150 deg., 0.5mm tip

NO IMAGE  
AVAILABLE

**P12839 : OCHSNER TROCAR, 14 FR**

For 10 Fr catheter