

## Hand-Held Retractors, Hooks

**054135 : COTTLE SKIN HOOK**

Single prong, 4-7/8" (12.5 cm)

**KM48217 : Inge Retractor**

with teeth, 6 1/2" (16.5cm)

**054110 : JOSEPH SKIN HOOK**

Sharp, 5.6 mm hook, double-pronged (10 mm between prongs), 6-1/8" (15.6 cm)

**054600 : Ragnell Retractor**

Double ended, Blade A 8 x 3mm / Blade B 15 x 5mm, 5 3/4" (14.5cm)

**054142 : KLEINERT-KUTZ SKIN HOOK, MEDIUM**

5 mm hook, single prong, 6-1/4" (16 cm)

**KM48218 : Inge Retractor**

without teeth, 6 1/2" (16.5cm)



**351784 : Henly Deep Center Blades Retractor**

5/8" (1.6cm) width x 1 1/2" (3.8cm) / 1 3/4" (4.5cm) depth blade, 4 7/8" (12.5cm)




**442104 : JOSEPH SKIN HOOK**

5.5 mm hook, double-pronged (5.5 mm between prongs), 6-1/4" (15.8 cm)



**164754 : Children's Retractor**

Blades 20mm long x 9.5mm wide (3/4 x 3/8) / 20mm long x 6mm wide (3/4 x 1/4), 6 1/8" (15.5cm)



**355175 : COARCTATION HOOK**

Blunt, right-angled, 7 mm hook, 10-1/8" (25.7 cm)



**442295 : Freeman Face Lift Retractor**

1 5/8" (4.0cm) width blade, 7 1/8" (18.0cm)



**KM52055 : KMEDIC SKIN HOOK**

Sharp, 2 mm hook, 4-3/4" (12.1 cm)




**060120 : SISSON NERVE HOOK**

Blunt, 7 mm hook, 90 angled, 6-1/8" (15.5 cm)



**KM48230 : Inge Retractor**

with teeth, 10 3/8" (26.5cm)



**243201 : Lowsley Prostatic Tractor**

Straight, 12 7/8" (33.0cm)



**654181 : ADSON DISSECTOR HOOK**

Sharp, 8" (20.3 cm)



**654191 : Frazier Dura Hook**

654191: Sharp, 4-7/8" (12.3 cm)



**KM46848 : Hohmann Retractor**

22mm width blade, 3-holes, 11" (28cm)




**KM46907 : Fukuda Style Shoulder Retractor, Small**

32mm width x 81mm length blade, 7" (18.0cm)



**KM48343 : BONE HOOK, MEDIUM**

Sharp, 11 mm x 18 mm hook, 9" (22.9 cm)



**341635 : Brewster Phrenic Retractor Nested Pair**

341635 Nested set of 2 Retractors: Outside Retractor (double-ended) - overall length 149 mm (5-7/8") -Large End - 13 mm wide x 35 mm deep -Small End - 10 mm wide x 24 mm deep Inside Retractor (double-ended) - overall...



**KM48332 : BONE HOOK**

Large, sharp, 28 mm hook, 5.0 mm diam., 7-7/8" (20 cm)




**164750 : Senn Retractor**

Sharp, 3-pronged, 1/4" (6mm) width x 3/4" (2.0cm) depth blade, 6 1/4" (16.0cm)



**841400 : K-P HOOK**

Blunt, 13" (33 cm)



**054120 : COTTLE SKIN HOOKS**

Sharp, 7 mm hook, double-pronged (12 mm between prongs), 5-7/8" (14.8 cm)