


Hand-Held Retractors, Hooks**341773 : Harrington Splanchnic Retractor**

Hollow handle, 2 1/2" (6.5cm) width x 5" (12.7cm) depth blade, 12" (30.5cm)

**164715 : Army-Navy Retractor**

5/8" (1.5cm) width blade, 8 1/2" (21.5cm), 2/set

**442145 : Gillies Skin Hook, Small, Delicate**

442145: 2.1 mm hook diameter (small). 0.7 mm hook width (delicate).

**164540 : Richardson Appendectomy Retractor**

3/4" (1.9cm) width x 2" (5.0cm) depth blade, 9 1/2" (24.0cm)

**442106 : JOSEPH SKIN HOOK**

7 mm hook, double-pronged (7 mm between prongs), 6-1/4" (15.9 cm)

**KM46827 : Hohmann Retractor**

15mm width blade, 6" (15.0cm)



164752 : Senn Retractor

Blunt, 3-pronged, 1/4" (6mm) width x 3/4" (2.0cm) depth blade, 6 1/4" (16.0cm)




KM69790 : JOHNSON-TUCKER HOOK AND FORK

Double-ended, delicate, 6-1/2" (16.5 cm)



351790 : Cooley Mitral Valve Retractor, Large

Curved right, 1 7/8" (4.7cm) width x 1 3/8" (3.5cm) depth blade, 7 3/8" (19.0cm)



KM52040 : FREER SKIN HOOK

Sharp, 2.5 mm hook, double-pronged, 6" (15.2 cm)




164650 : Brewster Retractor

3/4" (2.0cm) width x 1 3/4" (4.5cm) depth blade, 8 1/2" (21.5cm)



357532 : MICRO VASCULAR NERVE HOOK

Small, 7-1/4" (18.4 cm)



054130 : GILLIES CONVERSE SKIN HOOK, LARGE

Sharp, 4 mm hook, 7-1/8" (18 cm)



351788 : Cooley Mitral Valve Retractor, Small

Curved right, 7/8" (2.2cm) width blade, 7 1/4" (18.5cm)




728175 : CURVED J-HOOK SUCTION COAGULATOR ELECTRODE

1-1/4" (3.3 cm)



KM69990 : Cox Metatarsal Spreader

11mm width blade, 4 1/2" (11.5cm)



341780 : Wylie Splanchnic Retractor

Hollow handle, 2 1/2" (6.3cm) width x 6 3/4" (17.1cm) depth blade



422250 : O'CONNOR HOOK

4-7/8" (12.3 cm)




054165 : TRACHEA HOOK

Sharp, 8.7 x 1.4 mm hook, 7-1/2" (19 cm)



212295 : Pryor-Pean Retractor

1" (2.5cm) width x 4 1/2" (11.5cm) depth blade, 10 3/8" (26.5cm)



054107 : JOSEPH SKIN HOOK

Sharp, 5.6mm hook, Double-pronged (6mm between prongs), 6 1/4" (15.9cm)



351789 : Cooley Mitral Valve Retractor, Medium

Curved right, 1 3/8" (3.4cm) width blade, 7 1/2" (19.0cm)



594402 : HOEN HOOK

Long prong, 9-1/4" (23.5 cm)



728150 : Hasson "S" Retractor

10 mm wide, malleable, Set of two



KM46906 : Fukuda Style Shoulder Retractor, Large

38mm width x 81mm length blade, 7" (18.0cm)