


## Table-Mounted Retractors




**168251 : Circuitract Round Ring, Medium, 10-1/2"**




**168252 : Circuitract Round Ring, Large, 12-1/2"**



**168253 : Circuitract Oval Ring, Small**



**168254 : Circuitract Oval Ring, Medium**



**168255 : Circuitract Oval Ring, Large**



**168256 : Circuitract Segmented Ring Oval, Small**



**168257 : Circuitract Segmented Ring Oval, Medium**



**168258 : Circuitract Segmented Ring Oval, Large**



**168259 : Circuitract Ring Segments, 4"**



**168260 : Circuitract Ring Segments, 6"**



**168261 : Circuitract Ring Segments, 8"**



**168262 : Circuitract Bariatric Segmented Ring, Extra-Large**



**168263 : Circuitract Bariatric Vertical Post, 18-3/8"**




### **341100 : Strong T-Bar Retractor System**

Single post, table-mounted IMA Retractor with Cheanvechai Adjustable Rakes. (Shown with optional Retractor Light Beam, 163701). Consists of the following components:  
341112: Strong Support Bar 341113: Strong T-Bar 341114:...



### **341112 : Strong T Support Bar**



### **341113 : Strong T-Bar**




### **341114 : Pilling® Poly-Clamp Rail Clamp**

The Pilling Poly-Clamp is a highly versatile operating table rail clamp for use with the Strong T-Bar, Gomez, Poly-Tract Fowler and other retractor systems, as well as with other equipment to be attached to the operating...



### **341123 : Ochner-Favaloro Captive Washer**



### **341124 : Ochner-Favaloro Adjustable Rake & Arm**

Curved, 4 prongs



### **341126 : Ochner-Favaloro Knob & Clamp**



**341127 : Ochner-Favaloro Adjustable Knob**



**341129 : Ochner-Favaloro Micro Rake and Clamp**

Rake Assembly for Ochser-Favaloro Self-Retaining Retractors 341120 and 341128  
Curved, 4 prongs, 7/8" wide (2 cm)



**341131 : Cheanvechai-Favaloro Rake Assembly**

Curved, 4 prongs



**341133 : Cheanvechai-Favaloro Swivel Rake Assembly**

Curved, 4 prongs



**P14213 : Diesel Retractor Clamp**

