

Cannulas



NO IMAGE
AVAILABLE

P25394 : OSHER "J" CANNULA

27 Ga.



NO IMAGE
AVAILABLE

P26412 : SHORT, STRONG NEEDLE

Curved, 6 1/2" (16.5cm)



NO IMAGE
AVAILABLE

P26413 : LIGHT NEEDLE

Curved, 6 1/2" (16.5cm)



NO IMAGE
AVAILABLE

P26444 : BAYONET NEEDLE

Fine, 6 1/2" (16.5cm)



NO IMAGE
AVAILABLE

P26781 : HOUSE BARON SUCTION TUBE

Extra long, 5 Fr., 110mm



NO IMAGE
AVAILABLE

P26782 : HOUSE BARON SUCTION TUBE

Extra long, 7 Fr., 110mm

NO IMAGE
AVAILABLE

P26786 : BRACKMAN SUCTION TUBE

5 Fr., 3" (7.5cm) working length

NO IMAGE
AVAILABLE

P26788 : BRACKMAN SUCTION TUBE

7 Fr., 3" (7.5cm) working length

NO IMAGE
AVAILABLE

P26795 : BELLUCCI EAR SUCTION TUBE

Tube and 22 Gauge tip only

NO IMAGE
AVAILABLE

P26796 : BELLUCCI EAR SUCTION TUBE

Tube and 24 Gauge tip only

NO IMAGE
AVAILABLE

P26812 : ROSEN SUCTION TUBE

14 Gauge

NO IMAGE
AVAILABLE

P26813 : ROSEN SUCTION TUBE

15 Gauge

NO IMAGE
AVAILABLE

P26814 : ROSEN SUCTION TUBE

16 Gauge




P26823 : SCHUKNECHT SUCTION TUBE

20 Gauge



P26841 : HOUSE SUCTION IRRIGATOR

2.5 x 4 Fr., 95mm



P26843 : HOUSE SUCTION IRRIGATOR

4 x 5 Fr., 95mm



P26844 : HOUSE SUCTION IRRIGATOR

5 x 7 Fr., 95mm




P26847 : HOUSE SUCTION IRRIGATOR

8 x 12 Fr., 95mm



P26848 : HOUSE SUCTION IRRIGATOR

12 x 14 Fr., 95mm



P27806 : SINUS SUCTION TUBE

Lightly curved, 3.4mm, 6 1/8" (15.5cm)

AVAILABLE

NO IMAGE
AVAILABLE

P27817 : SINUS SUCTION TUBE, LONG

2.5mm

NO IMAGE
AVAILABLE

P27818 : SINUS SUCTION TUBE, SHORT

2.5mm

NO IMAGE
AVAILABLE

P27820 : TROCAR AND CANNULA

4mm, fenestrated

NO IMAGE
AVAILABLE

P27821 : TROCAR AND CANNULA

5mm, oblique

NO IMAGE
AVAILABLE

P27831 : SEPTUM NEEDLE

22 Ga.