


**P17963 : MICRO DIAMOND DUST NEEDLE HOLDER**


Lightweight, fine, w/lock, 7 1/4" (18.5cm)

**351774 : Crawford Aortic Retractor**


22mm width x 3/4" (2.0cm) depth blade, 7 7/8" (20.0cm)

**355098 : Garrett Vascular Dilator**Long 4.5 mm, malleable shaft *Related: 355080 Short 1.0 mm*
355081 Short 1.5 mm◆ Long 1.0 mm
◆ Long 1.5 mm 355082 Short 2.0...**152850 : WANGENSTEEN NEEDLE HOLDER**

16.5 x 2.3mm jaw, 10 3/4" (27.5cm), TC

**355311 : COOLEY DILATOR**

1.5mm, 5" (12.5cm)

**P17968 : MICRO DIAMOND DUST NEEDLE HOLDER**

w/ lock, 8 1/4" (21cm)



351801 : DEBAKEY TISSUE FORCEPS

Straight, Delicate, 1mm tip, 6" (16 cm)




355096 : Garrett Vascular Dilator

Long 3.5 mm, malleable shaft *Related: 355080 Short 1.0 mm*
355081 Short 1.5 mm


◆ Long 1.0 mm

◆ Long 1.5 mm *355082 Short 2.0...*




359020 : Potts Micro Titanium Scissors

Angled 130, 6 3/4" (17.0cm)



354985 : PILLING INTRACARDIAC NEEDLE HOLDER


17 x 1.8mm jaw, 7" (18cm), TC



355306 : DeBakey® Vascular Dilator

10.0 mm, 7-1/2" (19 cm) long *Related: 355295 1.0 mm, 7-1/2" (19 cm) long*
Set of 12, 15" (38 cm), in case 355296 1.5 mm, 7-1/2" (19 cm) long...

355290



359017 : Potts Micro Titanium Scissors

Angled 45, 7 1/4" (18.5cm)




841051 : KAISER-PILLING ATRAUGRIP RING FORCEPS (5AG)

Reverse curved, DeBakey type jaw, 11-3/4" (30 cm)



342585 : PILLING FORCEPS

Slight Curved, Serrated, 2" (5.1cm) tip length, 2mm tip width, 9 1/4" (23.5cm)



342721 : WILLAUER ALLIS TISSUE FORCEPS

Straight, 5 x 6 Teeth, 11 1/4" (29 cm)



354973 : Castroviejo Needle Holder

Round handle and lock. Diama-grip jaws. Curved, jaw approximately 11.5 x 0.9 mm, 6" (15 cm)




355097 : Garrett Vascular Dilator

Long 4.0 mm, malleable shaft *Related: 355080 Short 1.0 mm*  Long 1.0 mm
355081 Short 1.5 mm  Long 1.5 mm *355082 Short 2.0...*



351752 : Blanco Valve Retractor

7mm width x 3/8" (1.0cm) depth blade, 11" (28.0cm)



355300 : Debakey® Vascular Dilator

4.0 mm, 7-1/2" (19 cm) long *Related: 355295 1.0 mm, 7-1/2" (19 cm) long* *355290*
Set of 12, 15" (38 cm), in case 355296 1.5 mm, 7-1/2" (19 cm) long...



341531 : UNIVERSAL CLAMP FOR WILCOX PEDIATRIC RETRACTOR

Set of 2




351750 : Blanco Valve Retractor

5mm width x 3/8" (1.0cm) depth blade, 11 1/2" (29.0cm)




352132P : Beall Micro Coronary Scissors

Angled 25, 7" (18.0cm)




341514 : Chaux-Lemole Bi-Angle T-Frame for Cooley and Chaux Retractors

Bi-Angle T-Frame, for additional third point atrial retraction during mitral valve surgery, fits Cooley-style and Chaux Retractors. Arm tip maximum diameter 0.252", overall width 9-7/8" (25 cm) *Related: 341150 Chaux...*



354986 : PILLING INTRACARDIAC NEEDLE HOLDER

14.5 x 1.6mm jaw, 8" (20.5cm), TC



355085 : Garrett Vascular Dilator

Short 3.5 mm, malleable shaft *Related: 355080 Short 1.0 mm* ◆ *Long 1.0 mm*
355081 Short 1.5 mm ◆ *Long 1.5 mm* *355082 Short 2.0...*

