


**506737 : DREW ALLIGATOR CLIP APPLYING FORCEPS**


Straight, serrated, 2 mm tip, 9-7/8" (25 cm)

**341740 : DeBakey®-Cooley Retractor**

Hollow handle, 2 1/4" (5.7cm) width x 7 1/4" (18.4cm) depth blade, 12 1/4" (31.0cm)

**342594 : LEJEUNE FORCEPS**

Slight Curved, Serrated, 1 3/8" (3.5cm) tip length, 1.8mm tip width, 9 3/8" (24cm)

**342655 : DEBAKEY LUNG GRASPING FORCEPS**

Straight, 9" (23 cm)

**357645 : MICRO RING TIP FORCEPS**

Smooth, 1.5mm tip, 7" (18 cm)

**341147 : CHAUX BLADE ASSEMBLY & VALVE RETRACTOR**



**341678 : Allison Lung Retractor, Child**

1 1/2" (4.0cm) width blade, 10 1/4" (26.0cm)



**351748 : CARDIAC VALVE FORCEPS**

Serrated, 3mm tip, 9 1/2" (25 cm)



**341148 : CHAUX RETRACTOR BLADE ASSEMBLY**

Side Blade 1 3/4" (4.5cm) deep, Extra-long



**341154 : Chaux Thumb Screw for Blade Assembly**



**341155 : Chaux Sternal IMA-Mitral Valve Retractor Set**

Three retractors in one, can be used as a sternal/IMA and valve retractor. Includes four sternal blades with individual and independent dual movement. Sliding sternal blade assemblies allow movement on horizontal planes of...



**341156 : Chaux Take-Up for Retractor Arm**



**352298 : Haimovici Scissors**

Acutely Curved, 10" (25.5cm), Delicate




**355100 : Garrett Dilator Holder**

8-1/8" (20.5 cm) *Related: 355080 Short 1.0 mm* ♦ Long 1.0 mm 355081 Short 1.5 mm ♦ Long 1.5 mm 355082 Short 2.0 mm...



**342160 : DeBakey-Metzenbaum Scissors**

Curved, 5 3/4" (14.5cm), TC




**355305 : DeBakey® Vascular Dilator**

9.0 mm, 7-1/2" (19 cm) long *Related: 355295 1.0 mm, 7-1/2" (19 cm) long 355290 Set of 12, 15" (38 cm), in case 355296 1.5 mm, 7-1/2" (19 cm) long...*



**341177 : DeBakey Chest Retractor, Adolescent w/2 Sets of Aluminum Side Blades**

Spread 5 1/2" (14cm), Side Blades 1 3/16" (3.0cm) wide x 1 1/4" (3.2cm) deep & 1 3/16" (3.0cm) wide x 1 3/4" (44cm) deep



**354996 : RYDER EXTRA FINE NEEDLE HOLDER (FRENCH EYE)**

14.5 x 1.2mm jaw, 8" (20.5cm), TC



**341198 : Bonchek-Carpentier Atrial Retractor**

Spread 6 7/8" (17.5cm), Side Blade 3 15/16" (10cm) wide x 1 1/4" (3.2cm) deep




**341205 : Replogle Retractor, Infant**

Spread 2 15/16" (7.5cm), Side Blade 9/16" (1.5cm) wide x 5/8" (1.6cm) deep, Arm length 5 7/8" (15cm)



### 342569 : DEBAKEY FORCEPS

Right Angle, Longitudinal Serrations, Cross-serrated tip, 7/8" (2.15cm) tip length, 11 3/8" (29cm)



### 341224 : Farrell Rubio Retractor W/Interchangeable Blades

Adult, 2 sets of blades Blades 1-5/8" deep x 2-5/8" wide (4.1 x 6.7 cm) and 2-1/2" deep x 2-5/8" wide (6.3 x 6.7 cm) Spread 7-1/2" (19 cm) Arms curved 50 degrees



### 341305 : Amato Neonatal Sternal Retractor, large

Specially designed infant sternal retractors to reduce post-operative sternal deformity by dispensing pressure over greater surface area. Selection of 3 sizes available for proper fit. 6 blade, large Spread 2-1/2" (6.5 cm)...



### 341387 : TRANSPLANT DONOR STERNAL RETRACTOR, INFANT/CHILD

Spread 4 3/8" (11cm), Sharp, Side Blade 1/2" (1.2cm) wide x 1/2" (1.2cm) deep



### 341502 : LEMOLE-PILLING® VALVE RETRACTOR - Set #2 Fits Morse-style Retractors

Provides self-retaining valve retraction. Fits directly onto arms of many types of sternal retractors. The T-Frame cross bar attaches to the sternal retractor by means of a "C" clamp. 8" (20.3 cm) T-Frame cross bar allows...

