



341503 : LEMOLE-PILLING® VALVE RETRACTOR - Set #3 Fits Ankeney Retractor

Provides self-retaining valve retraction. Fits directly onto arms of many types of sternal retractors. The T-Frame cross bar attaches to the sternal retractor by means of a "C" clamp. 8" (20.3 cm) T-Frame cross bar allows...




351814 : DEBAKEY TISSUE FORCEPS

Angled, 2mm tip, 9 1/2" (24 cm)




341515 : Lemole Side Arm Extension For Cooley And Chaux Retractors

For third point atrial retraction during mitral valve surgery, fits Cooley-style and Chaux Retractors. *Related: 341150 Chaux Multipurpose Sternal IMA Retractor 341155 Chaux Sternal/IMA Valve Retractor Set 341194,...*



351755 : Blanco Valve Retractor

9mm width x 1/2" (1.0cm) depth blade, 11 1/2" (29.0cm)




342606 : MIXER FORCEPS

Strongly Curved, Serrated (3cm) tip length, 10 3/4" (27.5cm)




152827 : MICRO SURGERY NEEDLE HOLDER

Curved, Flat handle, With lock, 12.8 x 1.4mm smooth jaw, 6" (14.5cm), TC




342166 : DeBakey-Metzenbaum Scissors

Curved, 9 1/8" (23.0cm), TC




353040 : SCHWARTZ CLAMP HOLDING FORCEPS

Angled, 7 1/2" (19 cm)




342725 : BABCOCK TISSUE FORCEPS

Serrated, 9 1/2" (24 cm)



352295 : Satinsky Scissors

Curved, 9 3/4" (25.0cm)




506748 : MEDIASTINAL BLUNT TIPPED DISSECTING FORCEPS

Serrated, 16 mm tip, 12" (31 cm), working length 7-1/2" (19 cm)



359048 : GERALD TITANIUM TISSUE FORCEPS

Serrated, 6 7/8" (18 cm)



343027 : STRATTE NEEDLE HOLDER

Reverse curved, 9" (23cm), TC

NO IMAGE
AVAILABLE

354908 : MATCHED SET OF FOUR DEBAKEY NEEDLE HOLDERS

Consists of 354900, 354902, 354904, 354906

357675 : Micro Vascular Round Handle Needle Holder

Smooth, stainless steel jaws, round handles. Curved, .6 mm jaws, with lock, 9" (23 cm)

355101 : Garrett Dilator Holder

5-1/8" (13 cm) *Related: 355080 Short 1.0 mm*

◆ Long 1.0 mm 355081 Short 1.5 mm

◆ Long 1.5 mm 355082 Short 2.0 mm...

357677 : Micro Vascular Round Handle Needle Holder

Smooth, stainless steel jaws, round handles. Straight, 1 mm tungsten carbide jaws, 8" (20.5 cm)

342560 : DEBAKEY FORCEPS

C Curved, 3/8" (1cm) tip length, 8 1/8" (20.5cm)

355309 : COOLEY DILATOR

0.5mm, 5" (12.5cm)

359058 : MICRO TITANIUM NEEDLE HOLDER

Curved, Flat handle, With lock, 1mm Diama-Grip jaw, 7 1/8" (18cm)




359059 : MICRO TITANIUM NEEDLE HOLDER

Straight, Flat handle, With lock, 1.1mm Diama-Grip jaw, 8 1/4" (21cm)



506742 : SCOVILLE ALLIGATOR CLIP APPLYING FORCEPS

Straight, 2 mm tip, 9-7/8" (25 cm)



359061 : MICRO TITANIUM NEEDLE HOLDER

Straight, Flat handle, With lock, 0.9mm Diama-Grip jaw, 6 1/4" (16cm)



359062 : MICRO TITANIUM NEEDLE HOLDER

Curved, Flat handle, With lock, 0.9mm Diama-Grip jaw, 6 1/4" (16cm)



351907 : GERALD TISSUE FORCEPS

1 x 2 Teeth, 9 3/4" (25 cm)

