

**342006 : STERNAL LOCKING PLATES**

A device used to fixate size 5 (B&S Gage 20) monofilament stainless steel suture wire to achieve sternal approximation. Only 316L stainless steel sternal closure wire should be used with the locking plate. Suture wire of...

**342015 : GRUNWALD SUTURE RING**

13 1/4" (33.6cm)

**342020 : LEBSCHE STERNUM KNIFE**

Straight tip, Straight shaft, 9 7/8" (25cm)

**342040 : SWEET MODIFIED STERNAL PUNCH**


10" (25.5cm)

**342043 : Kirklin Sternal Awl**


Working length 2" (5.0cm), Pointed needle end, Wrench included, 7 1/8" (18.2cm)

**342044 : KIRKLIN NEEDLE ONLY FOR 342043**


2-3/4" (6.8 cm)

**342102 : SAUERBRUCH RONGEUR**

Double action. Square bite, straight jaws, 11-3/4" (30 cm). 13 mm *Related: 342105P 17 mm 342107 20 mm*

**342105P : SAUERBRUCH RONGEUR**


Double action. Square bite, straight jaws, 11-3/4" (30 cm). 17 mm *Related: 342102 13 mm 342107 20 mm*

**342107 : SAUERBRUCH RONGEUR**

Double action. Square bite, straight jaws, 11-3/4" (30 cm). 20 mm *Related: 342102 13 mm 342105P 17 mm*

**342141 : URSCHEL FIRST RIB RONGEUR**


Double action, Straight, 5 x 12mm, Working length 9" (23cm)

**342142 : DALE FIRST RIB RONGEUR**


Double action, Curved, 3 x 15mm, 14 1/2" (37cm)

**342143 : URSCHEL-LEKSELL RONGEUR**

Double action, Curved, 8.5 x 15mm, 12 1/4" (31cm)


**342501 : RUMEL FORCEPS**

A Curved, DeBakey Jaw, 2 3/4" (7.1cm) tip length, 9 5/8" (24.5cm)



342561 : DEBAKEY FORCEPS

D Curved, 1/2" (1.4cm) tip length, 2.8mm tip width, 8 1/4" (21cm)



342625 : RIENHOFF SWAN NECK CLAMP

Serrated, Curved shaft, 10 5/8" (27cm)




342632 : HARRINGTON CLAMP

Heavy construction Right Angle, Longitudinal Serrations, 1 1/2" (3.85cm) tip length, 12 3/8" (31.5cm)




342633 : HARRINGTON CLAMP

Light construction Right Angle, Longitudinal Serrations, 1 1/2" (3.85cm) tip length, 12 3/8" (31.5cm)



342710 : SEMB DISSECTING FORCEPS

Lightly Curved, Smooth, 10 1/4" (26 cm)



342800 : LEES BRONCHUS CLAMP

Right Angle, DeBakey Jaw, 1 3/4" (4.5cm) tip length, 8 1/8" (20.5cm)




342802 : LEES BRONCHUS CLAMP

Angled, DeBakey Jaw, 1 5/8" (4cm) tip length, Curved shaft, 10" (25.5cm)



342804 : BRONCHUS CLAMP

Angled, DeBakey Jaw, 3 1/4" (8.3cm) tip length, 9" (23cm)



342822 : SAROT BRONCHUS CLAMP

8.9cm tip length, 9 1/4" (23.5cm)



342823 : SAROT BRONCHUS CLAMP

3 1/2" 8.9cm tip length, 9 1/4" (23.5cm)



342950 : NELSON TISSUE FORCEPS

6 x 7 Teeth, 9" (23 cm)



343145 : ALEXANDER PERIOSTEOTOME

Double ended, Stainless steel blades. Cross-serrated handle, 8 1/4" (21cm)

