



## CVT

**343148 : DOYEN RASP**

Infant size, Single ended, Right side, 5 7/8" (15cm)

**343149 : DOYEN RASP**

Infant size, Single ended, Left side, 5 7/8" (15cm)

**343152 : DOYEN RASP**


Adult size, Single ended, Right side, 6 3/4" (17cm)

**343153 : DOYEN RASP**


Adult size, Single ended, Left side, 6 1/2" (16.5cm)

**343590 : CHEST DRAIN PULLING CLAMP**

Straight, Serrated, 2" (5cm) tip length, 2mm tip width, 8 5/8" (22cm)


**351316 : FITCH OBTURATOR 26FR**

17 7/8" (45.5cm)



**351318 : FITCH OBTURATOR 30FR**

16 1/8" (41.0cm)



**351348 : PILLING WIRE TWISTER FORCEPS**

3.2mm, 5 3/8" (13.5cm)



**351349 : RUBIO WIRE TWISTER**

4.6mm, 5" (12.5cm), TC




**351351 : MORSE TUBING HOLDER FOR SIZE 3/16" AND 1/4" TUBING**

5 3/8" (13.5cm)



**351352 : MORSE TUBING HOLDER FOR 5/16" AND 3/8" TUBING**

5 3/8" (13.5cm)




**351353 : MORSE TUBING HOLDER FOR 7/16" AND 1/2" TUBING**

5 3/8" (13.5cm)



**351354 : MORSE TUBING HOLDER FOR 9/16" AND 5/8" TUBING**

5 3/8" (13.5cm)




**351355 : MORSE TUBING HOLDER FOR 11/16" AND 3/4" TUBING**

5 3/8" (13.5cm)



**351358 : PRESBYTERIAN HOSPITAL TUBING CLAMP**

Straight, 3/8" (1cm) tip length, 6 1/2" (16.5cm)




**351370 : NICHOLS TUBING HOLDER**

7 1/4" (18.5cm)




**351390 : GUNNAR ROLLER CLAMP**

Straight, 7 1/8" (18cm)




**351470 : JEHLE CORONARY CATHETER, STERILE**

Cuff size diameter 9/32" (7.1mm), Catheter diameter .045" (1.14mm) I.D., .1" (2.54mm) O.D., overall length 20 1/4" (51.5cm)



**351472 : JEHLE CORONARY CATHETER, STERILE**

Cuff size diameter 11/32" (8.7mm), Catheter diameter .045" (1.14mm) I.D., .1" (2.54mm) O.D., overall length 20 1/4" (51.5cm)



**351474 : JEHLE CORONARY CATHETER, STERILE**

Cuff size diameter 13/32" (10.3mm), Catheter diameter .068" (1.73mm) I.D., .136" (3.45mm) O.D., overall length 20 1/4" (51.5cm)

AVAILABLE

NO IMAGE  
AVAILABLE

**351476 : JEHLE CORONARY CATHETER, STERILE**

Cuff size diameter 15/32" (11.9mm), Catheter diameter .068" (1.73mm) I.D., .136" (3.45mm) O.D., overall length 20 1/4" (51.5cm)

NO IMAGE  
AVAILABLE

**351490 : JEHLE CORONARY CATHETER, 4-HEAD, STERILE**

Cuff sizes 9/32", 11/32", 13/32", 15/32". Four tubes (2 tubes .1" O.D. x .045" I.D., 2 tubes .136" O.D. x .068" I.D.), overall length 20 1/2" (52cm)



**351500 : DEBAKEY SUCTION TUBE**

9 mm tip, 11-3/8" (29.0 cm)



**351501 : SUCTION TUBE TIP FOR 351500**

3/8" (1.0 cm)



**351502 : DEBAKEY SUCTION TUBE TIP**

Curved, 3 mm tip, 12-1/4" (31.0 cm)