



CVT

**351744 : REVUELTA-GARCIA ALLIS CLAMP**

Half Curved, DeBakey Jaw, 10 7/8" (27.5cm)

**351747 : Revuelta Garcia Valve Leaflet Retractor**

Curved, 24.0mm width blade

**351783 : Lemole Atrial Retractor**

1 7/8" (5.0cm) width blade, 10 1/2" (26.5cm)

**351810 : DEBAKEY TISSUE FORCEPS**

Angled, 2mm tip, 6" (16 cm)

**351824 : DEBAKEY TISSUE FORCEPS**

Delicate, 9 1/2" (24 cm), TC

**351830 : GREGORY SUTURE STAY CLAMP**

Straight, 3 7/8" (10cm)



351831 : GREGORY SUTURE STAY CLAMP

Curved, 3 7/8" (10cm)



351835 : ADAMS TISSUE FORCEPS

Delicate, 1mm tip, 8" (21 cm), TC



351839 : CORONARY ULTRA FINE TISSUE FORCEPS

9"



351841 : SELMAN TISSUE FORCEPS

DeBakey type jaw, 8" (21 cm)



351843 : DEBAKEY BAYONET TISSUE FORCEPS

7 3/4" (20 cm)



351845 : DEBAKEY BAYONET TISSUE FORCEPS

Fine, 7 3/4" (20 cm)



351856 : POTTS SMITH DRESSING FORCEPS

Fine, Serrated, 10 1/4" (26 cm)



351872 : POTTS SMITH TISSUE FORCEPS

Fine, Serrated and 1 x 2 Teeth, 8 1/4" (21 cm)



351910 : GERALD ATRAUGRIP TISSUE FORCEPS

Angled, DeBakey type jaw, 6 7/8" (18 cm)



352598 : CRAWFORD SUTURE RING

10" (25.4cm)




352605 : VALVULOTOME ANTEGRADE CUT

5-7/8" (15 cm)




352610 : MILLS VALVULOTOME

Blade 2 mm long, 11-3/4" (30 cm)



352611 : In-Situ Venous Valvulotome Storage Case

A sturdy, stainless steel storage case that holds up to four valvulotomes. 14-5/8" (37 cm)



352612 : IN-SITU VENOUS VALVULOTOME

Blade 2 mm long, 9" (23 cm)



352613 : IN-SITU VENOUS VALVULOTOME

Blade 2 mm long, 11-3/4" (30 cm)



352614 : IN-SITU VENOUS VALVULOTOME

Blade 2 mm long, 9-7/8" (25 cm)



352615 : GERBODE VALVE DILATOR

1 3/4" (4.5cm) wide



352630 : TUBBS VALVE DILATOR

1 3/4" (4.5cm) wide



352680 : THE BUSH RETROGRADE VALVULOTOME

Set of 3 (2, 3, 4 mm heads, 1 each), 15-3/4" (40 cm)

