



Home // CVT

CVT

**352681 : THE BUSH RETROGRADE VALVULOTOME**

2 mm head, 15-3/4" (40 cm)

**352682 : THE BUSH RETROGRADE VALVULOTOME**


3 mm head, 15-3/4" (40 cm)

**352683 : THE BUSH RETROGRADE VALVULOTOME**


4 mm head, 15-3/4" (40 cm)

**352690 : THE BUSH RETROGRADE VALVULOTOME**

Set of 3 (2, 3, 4 mm heads, 1 each), 37-3/8" (95 cm)

**352691 : THE BUSH RETROGRADE VALVULOTOME**

2 mm head, , 37-3/8" (95 cm)

**352692 : THE BUSH RETROGRADE VALVULOTOME**

3 mm head, 37-3/8" (95 cm)



352693 : THE BUSH RETROGRADE VALVULOTOME

4 mm head, 37-3/8" (95 cm)



352705 : AMATO PEDIATRIC DILATORS, SET OF 10 PROBES W/CASE

During the repair of complex congenital heart disease, the surgeon is required in many instances to locate orifices, trace the course of abnormal vessels, and perform various maneuvers of dilation and measurement of...



352800 : BAILEY AORTIC VALVE RONGEUR

Straight, 3.5 x 7mm



352802 : BAILEY AORTIC VALVE RONGEUR

Bite down 25, 3.5 x 7mm, Working length 5 1/8" (13cm)



352804 : BAILEY AORTIC VALVE RONGEUR

Bite down 70, 3.5 x 7mm, Working length 5 1/8" (13cm)



352806 : BAILEY AORTIC VALVE RONGEUR

Bite down 90, 3.5 x 7mm, Working length 5 1/8" (13cm)



352808 : BAILEY AORTIC VALVE RONGEUR

Angled left 45, 3.5 x 7mm, Working length 5 1/8" (13cm)

352810 : BAILEY AORTIC VALVE RONGEUR

Angled right 45, 3.5 x 7mm, Working length 5 1/8" (13cm)



352825 : WYLIE CALCIFIED ATERIAL PUNCH

Curved, 8 1/4" (21cm), Punches 1.3mm Hole



353001 : JOHNS HOPKINS BULLDOG CLAMP

Straight, Serrated, 1 1/2" (4cm)



353002 : JOHNS HOPKINS BULLDOG CLAMP

Straight, Serrated, 1 7/8" (5cm)



353003 : JOHNS HOPKINS BULLDOG CLAMP

Straight, Serrated, 2 1/4" (5.5cm)



353011 : JOHNS HOPKINS BULLDOG CLAMP

Curved, Serrated, 1 1/2" (4cm)



353012 : JOHNS HOPKINS BULLDOG CLAMP

Curved, Serrated, 1 7/8" (5cm)





353014 : JOHNS HOPKINS BULLDOG CLAMP

Curved, Serrated, 2 1/2" (6.5cm)




353016 : JOHNS HOPKINS BULLDOG CLAMP

Curved, Serrated, 3" (7.5cm)



353025 : HESS BULLDOG CLAMP

Angled, Serrated, 1 1/8" (2.7cm) tip length, 2 1/8" (5.5cm)



353035 : BULLDOG CLAMP APPLYING FORCEPS

Angled, Smooth, 9 1/2" (24 cm)



353041 : SCHWARTZ MICRO BULLDOG CLAMP

Straight, 1" (2.5cm)