

CVT

**353042 : SCHWARTZ MICRO BULLDOG CLAMP**

Slight Curved, 1/4" (6.4mm) tip length, 1" (2.5cm)

**353043 : SCHWARTZ MICRO BULLDOG CLAMP**

Strongly Curved, 3/8" (10.2mm) tip length, 1" (2.5cm)

**353050 : DEBAKEY CROSS ACTION BULLDOG CLAMP**

Straight, 3/4" (2cm) tip length, 3 1/8" (8cm)

**353052 : DEBAKEY CROSS ACTION BULLDOG CLAMP**

Straight, 1 1/8" (3cm) tip length, 3 1/2" (9cm)

**353054 : DEBAKEY CROSS ACTION BULLDOG CLAMP**

Straight, 2" (5.1cm) tip length, 4 1/4" (11cm)

**353056 : DEBAKEY CROSS ACTION BULLDOG CLAMP**

Straight, 2 5/8" (6.8cm) tip length, 4 7/8" (12.5cm)



353060 : DEBAKEY CROSS ACTION BULLDOG CLAMP

Curved, 5/8" (1.6cm) tip length, 2 7/8" (7.5cm)



353062 : DEBAKEY CROSS ACTION BULLDOG CLAMP

Curved, 1" (2.5cm) tip length, 3 1/4" (8cm)



353064 : DEBAKEY CROSS ACTION BULLDOG CLAMP

Curved, 1 5/8" (4.1cm) tip length, 3 7/8" (10cm)



353066 : DEBAKEY CROSS ACTION BULLDOG CLAMP

Curved, 2 1/8" (5.3cm) tip length, 4 1/4" (11cm)




353070 : SANTULLI BULLDOG CLAMP

Angled, DeBaKey Jaw, 1 5/8" (4.1cm) tip length, 2 3/4" (7cm)



353075 : WYLIE LUMBAR BULLDOG CLAMP

Slight Curved, DeBaKey Jaw, 7/8" (2.1cm) tip length, 3" (7.5cm)




353077 : GREGORY CAROTID BULLDOG CLAMP

Angled, DeBaKey Jaw, 5/8" tip length, 4 1/2" (11.5cm)

353078 : GREGORY CAROTID BULLDOG CLAMP

Angled, DeBakey Jaw, 1/2" tip length, 4 1/2" (11.5cm)



353080 : CALMAN CAROTID BULLDOG CLAMP

Straight, DeBakey Jaw, 1" (2.4cm) tip length




353085 : CALMAN CAROTID BULLDOG CLAMP

Angled, DeBakey Jaw, 1 5/8" (4cm) tip length, 3 3/4" (9.5cm)



353086 : JAVID CAROTID SHUNT CLAMP

Curved, Longitudinal Serrations, 1 3/4" (4.45cm) tip length, 7 3/8" (18.5cm)



353087 : JAVID CAROTID SHUNT CLAMP

Curved, Longitudinal Serrations, 2 3/8" (6cm) tip length, 6 3/4" (17cm)



353100 : GLOVER BULLDOG CLAMP


Straight, DeBakey Jaw, 3/4" (1.9cm) tip length, 2 1/4" (5.5cm)



353102 : GLOVER BULLDOG CLAMP

Straight, DeBakey Jaw, 1 1/4" (3.3cm) tip length, 1/4" (0.5cm)





353105 : GLOVER BULLDOG CLAMP

Straight, DeBakey Jaw, 1 5/8" (4cm) tip length, 3 1/2" (9cm)



353110 : GLOVER BULLDOG CLAMP

Curved, DeBakey Jaw, 5/8" (1.6cm) tip length, 2" (5cm)



353112 : GLOVER BULLDOG CLAMP

Curved, DeBakey Jaw, 1" (2.6cm) tip length, 2 1/2" (6cm)



353115 : GLOVER BULLDOG CLAMP

Curved, DeBakey Jaw, 1 3/8" (3.4cm) tip length, 3 1/8" (8cm)



353120 : GLOVER BULLDOG CLAMP

Angled, DeBakey Jaw, 1 1/4" (3.3cm) tip length, 1 7/8" (5cm)