

CVT

**353123 : WOOD BULLDOG CLAMP**

Straight, Serrated, 3/8" (8mm) tip length, 1 3/4" (4.5cm)

**353124 : WOOD BULLDOG CLAMP**

Angled, Serrated, 3/8" (8mm) tip length, 1 5/8" (4cm)

**353125 : WOOD BULLDOG CLAMP**

Straight, Serrated, 1/2" (12mm) tip length, 1 7/8" (4.5cm)

**353126 : WOOD BULLDOG CLAMP**

Angled, Serrated, 1/2" (12mm) tip length, 1 3/4" (4.5cm)

**353127 : WOOD BULLDOG CLAMP**

Straight, Serrated, 3/4" (19mm) tip length, 2" (5cm)

**353128 : WOOD BULLDOG CLAMP**

Angled, Serrated, 3/4" (19mm) tip length, 1 7/8" (5cm)



353131 : JACOBSON BULLDOG CLAMP

Straight, Serrated, 1 7/8" (5cm)



353132 : JACOBSON BULLDOG CLAMP

Curved, Serrated, 1 5/8" (4.5cm)



353135 : SENNING FEATHERWEIGHT BULLDOG CLAMP - STRONG

Straight, Serrated, 1 3/4" (4.5cm)



353136 : SENNING FEATHERWEIGHT BULLDOG CLAMP MEDIUM STRONG

Straight, Serrated, 1 3/4" (4.5cm)



353137 : SENNING FEATHERWEIGHT BULLDOG CLAMP - MEDIUM LIGHT

Straight, Serrated, 1 3/4" (4.5cm)



353138 : SENNING FEATHERWEIGHT BULLDOG CLAMP - LIGHT

Straight, Serrated, 1 3/4" (4.5cm)



353140 : RENAL BULLDOG CLAMP 10MMX15MM

Straight

353142 : RENAL BULLDOG CLAMP 10MMX25MM

Straight




353150 : DIETHRICH BULLDOG CLAMP

Straight, Serrated, 3/8" (8mm) tip length, 1 3/4" (4.5cm)



353151 : DIETHRICH BULLDOG CLAMP

Straight, Serrated, 1/2" (12mm) tip length, 2" (5cm)



353152 : DIETHRICH BULLDOG CLAMP

Straight, Serrated, 3/4" (20mm) tip length, 2 3/8" (6cm)



353153 : DIETHRICH BULLDOG CLAMP

Angled, Serrated, 8mm tip length, 1 7/8" (5cm)



353154 : DIETHRICH BULLDOG CLAMP

Angled, Serrated, 12mm tip length, 1 7/8" (5cm)



353155 : DIETHRICH BULLDOG CLAMP

Angled, Serrated, 3/4" (20mm) tip length, 2 1/8" (5.5cm)





353182 : Debakey Ring Handle Bulldog Clamp

353182 A short, sturdy clamp intended for use in place of the longer sizes of "squeeze handle" bulldogs. It measures only 2-3/4" (7 cm) from the center point of the box lock to the end of the ring handles. Narrow DeBakey...



353183 : GREGORY RING HANDLE BULLDOG CLAMP

Straight, DeBakey Jaw, 1.6mm tip width, Straight shaft, 4 3/4" (12cm)



353186 : Debakey Ring Handle Bulldog Clamp

353186 A short, sturdy clamp intended for use in place of the longer sizes of "squeeze handle" bulldogs. It measures only 2-3/4" (7 cm) from the center point of the box lock to the end of the ring handles. Narrow DeBakey...



353187 : GREGORY RING HANDLE BULLDOG CLAMP

Half Curved, DeBakey Jaw, 1 1/8" (2.7cm) tip length, 1.5mm tip width, Curved shaft, 4 1/8" (10.5cm)



353189 : Debakey Ring Handle Bulldog Clamp

353189 A short, sturdy clamp intended for use in place of the longer sizes of "squeeze handle" bulldogs. It measures only 2-3/4" (7 cm) from the center point of the box lock to the end of the ring handles. Narrow DeBakey...