


**353190 : GREGORY RING HANDLE BULLDOG CLAMP**

Angled 45, DeBakey Jaw, 1 5/8" (4.2cm) tip length, 1.5mm tip width, Curved shaft, 4 3/8" (11cm)

**353192 : DeBakey Ring Handle Bulldog Clamp**


353192 A short, sturdy clamp intended for use in place of the longer sizes of "squeeze handle" bulldogs. It measures only 2-3/4" (7 cm) from the center point of the box lock to the end of the ring handles. Narrow DeBakey...

**353195 : DeBakey Ring Handle Bulldog Clamp**


353195 A short, sturdy clamp intended for use in place of the longer sizes of "squeeze handle" bulldogs. It measures only 2-3/4" (7 cm) from the center point of the box lock to the end of the ring handles. Narrow DeBakey...

**353203 : STYLET FOR RUMEL TOURNIQUET**

10 1/4" (26cm)

**353206 : STYLET FOR WECK STYLE RUMEL TOURNIQUET**

9 3/8" (24cm)

**353207 : STYLET CHILD SIZE RUMEL TOURNIQUET**

7" (18cm)




353490 : LELAND JONES CLAMP

Straight, DeBakey Jaw, 3 3/8" (8.7cm) tip length, 2.4mm tip width, 8 1/8" (20.5cm)



353491 : LELAND JONES CLAMP

Angled 15, DeBakey Jaw, 3 1/2" (8.9cm) tip length, 2.4mm tip width, 8 1/8" (20.5cm)




353492 : LELAND JONES CLAMP

Angled 30, DeBakey Jaw, 3 5/8" (9.05cm) tip length, 2.4mm tip width, 8" (20.5cm)



353494 : GREGORY PROFUNDA CLAMP

DeBakey Jaw, 2 1/4" (5.8cm) tip length, 1.7mm tip width, 4 3/4" (12cm)




353495 : GREGORY PROFUNDA CLAMP

DeBakey Jaw, 3 1/4" (8.3cm) tip length, 2.8mm tip width, 6 1/8" (15.5cm)



353496 : GREGORY PROFUNDA CLAMP

DeBakey Jaw, 4" (10.2cm) tip length, 2.8mm tip width, 7" (18cm)




353500 : DEBAKEY COARCTATION CLAMP

Straight, 3 1/8" (7.9cm) tip length, 2.8mm tip width, 9" (23cm)




353503 : DEBAKEY COARCTATION CLAMP

Straight, 3" (7.6cm) tip length, 2.8mm tip width, Angled shaft, 8 3/4" (22cm)




353510 : DEBAKEY PATENT DUCTUS CLAMP

Straight, 2 1/4" (5.7cm) tip length, 2.8mm tip width, Straight shaft, 8 1/4" (21cm)




353513 : DEBAKEY PATENT DUCTUS CLAMP

Straight, 2 3/8" (5.9cm) tip length, 2.8mm tip width, Angled shaft, 8 1/8" (20.5cm)




353515 : DEBAKEY RENAL ARTERY CLAMP

3 1/2" 9cm tip length, Curved shaft, 10 1/8" (25.5cm)



353522 : WIRTHLIN SPLENO RENAL SHUNT CLAMP

3 1/8" 7.9cm tip length, 1.7mm tip width, Curved shaft, 6 1/2" (16.5cm)




353524 : STONEY SPLENO RENAL SHUNT CLAMP

DeBakey Jaw, 3 1/4" (8.3cm) tip length, 2mm tip width, Curved shaft, 9 1/4" (23.5cm)




353528 : DALE FEM-POP CLAMP

DeBakey Jaw, 4 1/8" (10.4cm) tip length, 2.8mm tip width, Curved shaft, 7 3/8" (19cm)



353531 : CLEVELAND CLINIC RENAL CLAMP

Left Curved, DeBakey Jaw, 3 3/8" (8.6cm) tip length, 9 5/8" (24.5cm)



353532 : CLEVELAND CLINIC RENAL CLAMP

Right Curved, DeBakey Jaw, 3 3/8" (8.6cm) tip length, 9 5/8" (24.5cm)




353533 : CLEVELAND CLINIC BABY RENAL CLAMP

Left Curved, DeBakey Jaw, 2mm tip width, 5 7/8" (15cm)



353534 : CLEVELAND CLINIC BABY RENAL CLAMP

Right Curved, DeBakey Jaw, 2mm tip width, 5 7/8" (15cm)



353535 : DEBAKEY PERIPHERAL VASCULAR CLAMP

ANGLED, 3 3/4" (9.4cm) tip length, 2.8mm tip width, Curved shaft, 7 1/8" (18cm)