

**353536 : DIETHRICH CAVAL OCCLUDING CLAMP**


Curved, DeBakey Jaw, 4 1/4" (10.7cm) tip length, Curved shaft, 9 3/4" (25cm)

**353538 : HENLY SUBCLAVIAN ARTERY CLAMP**


4 1/8" 10.4cm tip length, 2.8mm tip width, 8 3/8" (21cm)

**353539 : CHERRY SUPRACELIAC CLAMP**


Curved, 2 x 3 AtrauGrip Jaw, Curved shaft, 12 1/4" (31cm)

**353540 : DEBAKEY "S" PV CLAMP**


Curved, 4 1/8" (10.47cm) tip length, 2.8mm tip width, Curved shaft, 7 7/8" (20cm)

**353541 : DALE FEM-POP CLAMP**

Curved, DeBakey Jaw, 3" (7.5cm) tip length, 2.4mm tip width, Curved shaft, 7" (17.5cm)


**353542 : DALE FEM-POP CLAMP**

Curved, DeBakey Jaw, 3 3/8" (8.5cm) tip length, 2.4mm tip width, Curved shaft, 7 1/8" (18cm)



353543 : DALE FEM-POP CLAMP

Curved, DeBakey Jaw, 3 1/2" (8.9cm) tip length, 2.4mm tip width, Curved shaft, 7 1/4" (18.5cm)




353545 : DEBAKEY AORTIC CLAMP

Curved, 3 7/8" (9.8cm) tip length, 2.4mm tip width, Curved shaft, 7 1/8" (18cm)



353546 : FAVALORO PROXIMAL ANASTOMOSIS CLAMP

Curved, DeBakey Jaw, 4 1/8" (10.4cm) tip length, 2.8mm tip width, 8" (20cm)



353550 : DEBAKEY PV CLAMP

Curved, 2.3mm tip width, 6 3/4" (17cm)




353560 : DEBAKEY BAINBRIDGE CLAMP

Straight, 2.4mm tip width, 7" (17.5cm)



353562 : DEBAKEY BAINBRIDGE CLAMP

Straight, 3 1/8" (8cm) tip length, 2.8mm tip width, 7 7/8" (20cm)



353565 : DEBAKEY BAINBRIDGE CLAMP

Curved, 2 3/8" (5.9cm) tip length, 2.4mm tip width, 7" (17.5cm)


353567 : DEBAKEY BAINBRIDGE CLAMP

Curved, 3 1/4" (8.1cm) tip length, 2.4mm tip width, 7 5/8" (19.5cm)




353568 : DEMOS TIBIAL ARTERY CLAMP

Straight, DeBakey Jaw, 2 3/8" (6cm) tip length, 2.1mm tip width, Angled shaft, 6 3/4" (17cm)




353569 : WYLIE CLOSER CAROTID CLAMP

DeBakey Jaw, 2 1/2" (6.3cm) tip length, 1.7mm tip width, Curved shaft, 6 1/8" (15.5cm)




353570 : WYLIE EXTERNAL CAROTID CLAMP

DeBakey Jaw, 3 3/4" (9.5cm) tip length, 2.8mm tip width, Curved shaft, 8" (20.5cm)



353571 : WYLIE INTERNAL CAROTID CLAMP

DeBakey Jaw, 2.7mm tip width, Curved shaft, 8 1/4" (21cm)



353580 : KAY ANASTOMOSIS AORTIC CLAMP

Curved, DeBakey Jaw, 1 1/4" (3.2cm) tip length, Angled shaft, 7 5/8" (19.5cm)



353581 : STRONG AORTA CLAMP

Curved, DeBakey Jaw, Straight shaft, 8" (20cm)






353582 : LEMOLE-STRONG AORTA CLAMP

Curved, DeBakey Jaw, Straight shaft, 7 7/8" (20cm)



353583 : STRONG AORTA CLAMP

Curved, DeBakey Jaw, Angled shaft, 8 1/8" (20.5cm)



353584 : LEMOLE-STRONG AORTA CLAMP

Curved, DeBakey Jaw, Angled shaft, 8" (20.5cm)



353585 : KAY LAMBERT AORTIC ANASTOMOSIS CLAMP

Curved, DeBakey Jaw, 3 7/8" (9.8cm) tip length, Angled shaft, 8 1/8" (20.5cm)



353586 : LAMBERT-KAY AORTIC ANASTOMOSIS CLAMP

Curved, DeBakey Jaw, 4 1/4" (10.8cm) tip length, Angled shaft, 8 1/8" (20.5cm)