

**353588 : MITCHELL AORTOMY CLAMP**


DeBakey Jaw, 3 1/8" (7.8cm) tip length, Curved shaft, 8" (20.5cm)

**353595 : ZANGER ABDOMINAL AORTA CLAMP**


Curved, DeBakey Jaw, 3 3/4" (9.6cm) tip length, Curved shaft, 8 1/8" (20.5cm)

**353596 : ZANGER ABDOMINAL AORTA CLAMP**


Curved, DeBakey Jaw, 4 1/8" (10.6cm) tip length, Curved shaft, 9 3/4" (25cm)

**353597 : ZANGER ABDOMINAL AORTA CLAMP**


Curved, DeBakey Jaw, 5 1/8" (13cm) tip length, Curved shaft, 11 1/4" (28.5cm)

**353600 : DEBAKEY ANEURYSM AORTA CLAMP**

Curved, 3 3/4" (9.4cm) tip length, Curved shaft, 10 3/8" (26.5cm)

**353602 : DEBAKEY ANEURYSM AORTA CLAMP**

Curved, 4" (10.2cm) tip length, Curved shaft, 10 7/8" (27.5cm)




**353605 : DEBAKEY ANEURYSM AORTA CLAMP**

Curved, 2 x 3 AtrauGrip Jaw, 4 7/8" (12.54cm) tip length, Curved shaft, 12 1/2" (32cm)




**353607 : DEBAKEY ANEURYSM AORTA CLAMP**

Curved, 2 x 3 AtrauGrip Jaw, 5 3/4" (14.6cm) tip length, Curved shaft, 13 1/4" (33.5cm)



**353609 : HOWARD-DEBAKEY AORTA CLAMP**

Curved, 4 1/8" (10.5cm) tip length, Curved shaft, 10 3/4" (27.5cm)



**353610 : HOWARD-DEBAKEY AORTA CLAMP**

Curved, 2 x 3 AtrauGrip Jaw, 4" (10.3cm) tip length, Curved shaft, 12 1/2" (32cm)




**353617 : MORRIS ASCENDING AORTA CLAMP**

2 x 3 AtrauGrip Jaw, 3 1/4" (8.4cm) tip length, Curved shaft, 8 1/8" (20.5cm)




**353618 : MORRIS ASCENDING AORTA CLAMP**

2 x 3 AtrauGrip Jaw, 3 3/4" (9.4cm) tip length, Curved shaft, 9 1/2" (24cm)




**353619 : BEALL-MORRIS ASCENDING AORTA CLAMP**

2 x 3 AtrauGrip Jaw, 4 1/2" (11.5cm) tip length, Curved shaft, 7 1/4" (18.5cm)



**353620 : HUFNAGEL AORTA CLAMP**

Curved, DeBakey Jaw, Angled shaft, 9 5/8" (24.5cm)




**353622 : HUFNAGEL AORTA CLAMP**

Straight, DeBakey Jaw, 4" (10.3cm) tip length, Angled shaft, 9 5/8" (24.5cm)



**353623 : BAHNSON AORTA CLAMP**

Curved, DeBakey Jaw, 4 1/8" (10.5cm) tip length, Angled shaft, 9" (23cm)




**353625 : DEBAKEY FULL CURVED AORTIC CLAMP**

2 1/8" 5.4cm tip length, Curved shaft, 9 1/2" (24cm)



**353626 : DEBAKEY FULL CURVED AORTIC CLAMP**

2 x 3 AtrauGrip Jaw, Curved shaft, 11 1/8" (28.5cm)



**353627 : DEBAKEY FULL CURVED AORTIC CLAMP**

2 x 3 AtrauGrip Jaw, 2 1/4" (5.72cm) tip length, Curved shaft, 12 1/2" (31.5cm)




**353630 : DEBAKEY CLAMP**

Acutely Curved, 1 1/4" (3.3cm) tip length, Curved shaft, 9 5/8" (24.5cm)




**353631 : DEBAKEY CLAMP**

Medium Curved, 1 1/4" (3.2cm) tip length, Curved shaft, 10 1/8" (25.5cm)



**353632 : DEBAKEY CLAMP**

Slight Curved, 1 1/4" (3.2cm) tip length, Curved shaft, 10 3/8" (26cm)




**353635 : DEBAKEY CLAMP**

Acutely Curved, 2 1/2" (6.35cm) tip length, Curved shaft, 10 3/8" (26cm)



**353639 : MITCHEL-ADAM CLAMP**

Curved, DeBakey Jaw, 1 3/8" (3.5cm) tip length, Curved shaft, 9 7/8" (25cm)



**353640 : DeBakey® Modified Semb Ligature Carrier**

The distal 1.5 mm of the jaws are a DeBakey-type Design to securely hold a ligature or umbilical tape for passage around the aorta. Jaws curve on 2.5 cm radius. Curved, DeBakey Teeth, 10-3/8" (26.5 cm)

